

## COLOR MATCHING SYSTEM MADE STRESS-FREE: ACCURATE COLOR MATCHING USING HB BODY'S SPECTROPHOTOMETER X RITE MA91 + OEM 99

HB BODY Using XRITES latest technology in colorimetry has made color matching easier than ever before. The spectrophotometer is a portable instrument which can take up to 1000 individual measurements that can be downloaded to process on the COLOR MATCHING SYSTEM SOFTWARE.

Is used in combination with the COLOR SOLUTION COLOR CHIPS and the COLOR SOLUTION software, for accurate color matching.

Following our color matching process, the results will be more accurate than ever.



Gently place the spectrophotometer on the panel to be measured.




Firmly put pressure on the spectrophotometer till all light indications turn green. (not too much pressure will turn the sensors orange, too much pressure will turn the sensors red)



You can either connect the spectrophotometer on the computer by USB cable or you can use the Bluetooth technology.




Upload the data taken to the CMS software in the particular order and rename if desirable.



Choose the best possible solution given. Find the color chip, on the COLOR SOLUTION COLOR CHIPS CABINET.



Confirm the matching with the panel



Select the formula from the COLOR SOLUTION SOFTWARE



Mix the formula accordingly.